Steri-Pak EU PLUS Bowie-Dick Test Pack, REF BD 113+

Directions For Use

INDICATIONS FOR USE

The Steri-Pak EU PLUS Bowie-Dick Test Pack is for daily monitoring of pre-vacuum steam sterilizers operating at 134° C. The BD 113+ model is packaged in an easy to open box and designed to operate in sterilizers with sub-atmospheric or trans-atmospheric cycles as defined in ISO 11140-4:2007.

COLOR CHANGE

The indicator sheet will change from blue to black when exposed to steam sterilization cycle with acceptable vacuum cycles.

CRITICAL PARAMETERS

Steam sterilization cycle functioning at 134° C for 3.5 – 4 minutes.

INSTRUCTIONS FOR USE

1. Run a shortened cycle (omitting the drying phase) to properly heat the sterilizer.
2. Place the test pack on the bottom of the sterilizer nearest the drain.
3. Run a Bowie-Dick test cycle as specified by the sterilizer manufacturer.
4. Run the test for 3.5 minutes at 134° C. The cycle time must never exceed 4 minutes.
5. A short dry time of 2 to 5 minutes will not affect the test.

INTERPRETATION

1. Remove the test pack from the sterilizer and open immediately. Pull easy open zip tab across box front to open.
2. Remove indicator sheet from the pack and dispose of remaining pack materials.
3. If the sterilizer is functioning normally, the indicator graphic will show a uniform color change (black) throughout the entire pattern.
4. When a significant amount of residual air is present, the indicator shows a blue area in the center.

SAFETY PRECAUTIONS

- *CAUTION: the test pack will be hot and should be opened carefully to avoid thermal injury.*

STORAGE

- Store at normal room temperature 50° - 100°F (10° – 38° C) and 10-70% Relative Humidity.
- Protect from moisture and excess humidity.
- Store exposed indicator sheet away from light sources to prevent reversion back to original blue color.

EXPIRY DATE

The expiry date is printed on the product packaging.

LOT NUMBER

A unique identification code, **LOT**, is printed on each indicator and packaging labels.

INTERFERING SUBSTANCES OR CONDITIONS

There are NO KNOWN INTERFERING SUBSTANCES OR CONDITIONS that could affect the intended use of the indicator or adversely affect the indicator performance.

RELEASE OF TOXIC SUBSTANCES

The indicator releases NO KNOWN TOXIC SUBSTANCES in sufficient quantities to cause either a health hazard, or a hazard to the intended properties of the product being sterilized before, during or after the sterilization process.

DECLARATION OF CONFORMITY

SteriTec Products Mfg. Co., Inc. declares that the Steri-Pak EU Plus Bowie Dick Test Pack Indicator conforms to all applicable portions ISO 11140-1:2005 and EN ISO-11140-1:2009 pertaining to Class 2 Indicators for use in specific tests.

Steri-Pak EU Plus Bowie Dick Test Pack have been tested during product validation according to EN ISO 11140-4:2007 and fully comply with the performance requirements after exposure to dry saturated steam when tested under the conditions specified in EN ISO 11140-4:2007.

This also certifies that the above testing was conducted in a steam sterilizer according to EN ISO 11140-4:2007 and using Air Injection Apparatus as described in EN ISO 11140-4:2007.

Manufactured By:
SteriTec Products Mfg. Co., Inc.
74 Inverness Drive E
Englewood, CO 80112 USA
(303) 660-4201 FAX: (303) 660-4213

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Fabbricante: SteriTec Products Mfg. Co., Inc., Englewood, CO 80112 (USA)