



# Product Technical Data Sheet

Form No.: FTP-009 Rev.: B

**Product Name:** SteriTec Dry Heat Indicator Label

**Product Type:** Process Indicator

**Reorder Number:** CI 113

**Description:** The SteriTec Dry Heat Indicator Labels are designed for use in a dry heat sterilizer operating at 160° C. The label provides a clear identification of processed versus unprocessed items. When used as directed, the Dry Heat Label gives visible confirmation that the item was exposed to a dry heat process.

**Performance Characteristics:**

When exposed to dry heat conditions, the indicator bar changes color.

**Color Change:** Green to Black.

**Complete Color Change Time:** 10 minutes +/- 5 minutes at 160° C.

**Indicator Stability:** Excellent before and after sterilization

**Physical Characteristics:**

**Dimensions:** 0.75" W X 1.14" L

**Indicator Substrate Material:** Removable adhesive label

**Quantity / Box:** 1000

**Boxes / Case:** 10

**Selling Unit:** Case

**Environmental Limitations:**

**Storage:** 50° - 100° F (10° - 38° C).

**Shelf Life:** 5 years from date of manufacture.

**FDA 510(k) Release to Market Date:** April 2, 1997 K953122

**Standards Classifications:**

AAMI ST60, Class 1

EN 867-1, Class A

ISO 11140-1, Class 1