



# Product Technical Data Sheet

Form No.: FTP-167 Rev.: B

**Product Name:** SteriTec Dry Checks Dry Heat Indicator Strip

**Product Type:** Process Indicator

**Reorder Number:** CI 136

**Description:** The SteriTec Dry Checks Dry Heat Indicator Strips are designed for use in a dry heat sterilizer operating at 160° C. The indicator provides a clear identification of processed versus unprocessed items. When used as directed, the Dry Checks Indicator gives visible confirmation that the item was exposed to a dry heat process.

**Performance Characteristics:**

When exposed to dry heat conditions, the indicator bar changes color.

**Color Change:** Green to Black.

**Complete Color Change Time:** 10 minutes +/- 2.5 minutes at 160° C.

**Indicator Stability:** Excellent before and after sterilization

**Physical Characteristics:**

**Dimensions:** 0.75" W X 3.8" L

**Indicator Substrate Material:** 80 lb Cover

**Quantity / Box:** 250

**Boxes/Case:** 20

**Selling Unit:** Case

**Environmental Limitations:**

**Storage:** 50° - 100° F (10° - 38° C).

**Shelf Life:** 3 years from date of manufacture.

**FDA 510(k) Release to Market Date:** April 2, 1997 K953122

**Standards Classifications:**

AAMI ST60, Class 1

EN 867-1, Class A

ISO 11140-1, Class 1