



Product Technical Data Sheet

Form No.: FTP-143-1 Rev.: A

Product Name: SteriTec International EZTest Incubator with Thermometer

Product Type: Incubator

SteriTec Reorder Number: IC 105

Customer Name: SteriTec

Customer Catalog Number: N/A

Description: The Biological Indicator Incubator is designed for single temperature incubation of either the EO and Steam EZTest Biological Indicators. The thermometer, included with the incubator, provides easy verification of incubation temperature.

Performance Characteristics:

Operating Range: Incubator: 37° – 57° C, Thermometer: 10° – 60° C (+/- .5° C).

Electrical Rating: Incubator: 50-60 Hz, 12 V, European Power Pack: 230 VAC, 50 Hz

Physical Characteristics:

The incubator has 14 wells for placement of the biological indicator vials and one well for placement of the thermometer. Switch on bottom panel for switching from 37°C to 57°C incubation temperatures.

Accessories: Thermometer (I14/T), Power Pack (P2-1000)

Materials: Incubator: Aluminium body, Thermometer: Glass Tube with alcohol.

Dimensions: Incubator: 3.5" W X 2" H X 4" D, Well Diameter = .375" Dia X 1" D
Thermometer: 3" L X .26" Dia

Units / Box: 1

Selling Unit: Box

Environmental Limitations:

Storage: Store at room temperature, 17° – 27° C. RH 20-80% non-condensing.

Shelf Life: N/A

SteriTec FDA 510(k) Release to Market Date: None. Industrial product. NOT a medical device.

Standards Classifications: N/A