



# Product Technical Data Sheet

Form No.: FTP-145 Rev.: A

**Product Name:** SteriTec Autoclave Tape Dispenser

**Product Type:** Tape Dispenser

**SteriTec Reorder Number:** TD 101

**Customer Name:** SteriTec

**Customer Catalog Number:** N/A

**Description:** The Autoclave Tape Dispenser was designed for dispensing 1" wide or less tape with a 3" core. For use with all types of Sterilization Tapes.

**Performance Characteristics:** N/A

**Physical Characteristics:**

**Materials:** Plastic

**Dimensions:** 7" L X 2" W X 3" H, 3" core, 1.25" tape trough.

**Units / Box:** 1

**Selling Unit:** Box

**Environmental Limitations:** N/A

**Shelf Life:** N/A

**SteriTec FDA 510(k) Release to Market Date:** None. Industrial product. NOT a medical device.

**Standards Classifications:** N/A