

# EU PLUS TEST PACK



## EU PLUS Bowie-Dick Test Pack with Blue-to-Black Ink Technology

Catalog No. **BD 113+**

- Superior ink chemistry aids interpretation  
Blue color remains if air leak is present
- Continuous ink design from center to edge  
Makes air leaks easier to confirm
- Permanent color change  
Can be used as a permanent record when stored  
away from light sources
- Lead and toxic heavy metal free
- Complete color change after 3.5 minutes at 134°C
- ANSI/AAMI/ISO 11140-5
- ISO 11140-1, class 2
- ISO 11140-4

### ORDERING INFORMATION

PART #	PACKAGED
BD 113+	1 CASE = 20 PACKS

**Steri-Tec**  
[www.steritecproducts.com](http://www.steritecproducts.com)