

BIOLOGICAL TEST PACK

STEAM



Disposable Steam Biological Test Pack with Instant Readout Integrator

Catalog No. BI 103/INT

- Test pack parallels performance of AAMI biological test packs

Disposable test pack provides reliable, consistent results test after test at a low cost

- Contains integrator card

Integrator card inserted in the test pack will not complete a purple-to-green color change until Sterility Assurance Level 10^{-6} (2x biological kill) has been met (3 minutes at 132°C-134°C)

Card provides a permanent record of the success of the cycle

Full cycle monitor that is sterilant specific meaning it must have steam to change

Card is lead and toxic heavy metal free

- Contains self-contained biological indicator
- Small, compact and easy to use with instant results
- Biological Indicator: ANSI/AAMI/ISO 11138
- Test Pack: ANSI/AAMI ST 79
- Integrator Card: FDA cleared integrator

ORDERING INFORMATION

PART #	PACKAGED
BI 103/INT	1 CASE = 30 PACKS AND 5 CONTROLS

SteriTec
www.steritecproducts.com