


CROSS-CHECKS EO



Inside-the-Pack/Pouch EO Integrator

to Check Time, Temperature, Humidity and EO Gas Concentration

Catalog No. CI 106

- Surpasses performance of biological spore strips

- Advanced accuracy

Requires all parameters for EO gas sterilization be met before color change occurs

Distinctive yellow-to-brown color change is not reversible

- Plastic lamination

Plastic lamination protects sterilized items from potential ink transfer

- Easy to interpret

Color reference arrow printed on each integrator

- FDA cleared integrator

- ISO 11140-1, class 4

ORDERING INFORMATION

PART #	PACKAGED
CI 106	1 CASE = 24 BAGS OF 250

SteriTec
www.steritecproducts.com