

CROSS-CHECKS



Inside-the-Pack Chemical Indicators for Pre-Vacuum Steam Sterilization

Catalog No. **CI 107, CI 111, CI 134**

- Surpasses performance of biological spore strips

Requires all parameters for steam sterilization to be met before color change can occur

Distinctive white-to-black (CI 107, CI 111) or blue-to-black (CI 134) color change is not reversible and indicator can be used as a permanent record

- Plastic lamination

Short version (CI 111) is plastic laminated to protect sterilized items from potential ink transfer

- ISO 11140-1, class 4

Checks proper time and temperature in a pre-vacuum steam sterilizer operating at 132°C (270°F) for US standards and 134°C (273°F) for international standards

- FDA cleared

ORDERING INFORMATION

PART #	PACKAGED
CI 107	1 CASE = 16 BAGS OF 250
CI 111	1 CASE = 24 BAGS OF 250
CI 134	1 CASE = 24 BAGS OF 250

SteriTec

www.steritecproducts.com