

STEAM CHARTS



Chemical Integrators 121°C/134°C

for Use in Testing Steam Sterilizers Operating at 121°C or 134°C

Catalog No. CI 108, CI 109, CI 108/121-134

- Full cycle monitors that are sterilant specific meaning it must have steam to change
- Temperature specific
 - CI 108 - testing at 121°C
 - CI 109 - testing at 134°C
 - CI 108/121-134 - testing for both 121°C or 134°C
- Performs equivalently to biological indicators at 121°C and 134°C
- Multiple color changing spots
 - First spot indicates if the item has been processed in a steam sterilizer
 - Second spot indicates if correct conditions for sterilization have been met
 - Central spot on all indicators shows excessive sterilization
- Lead and toxic heavy metal free
- Permanent color change
- ISO 11140-1, class 5

ORDERING INFORMATION	
PART #	PACKAGED
CI 108	1 CASE = 20 BAGS OF 250
CI 109	
CI 108/121-134	

SteriTec
www.steritecproducts.com