

EZESTEO



Biological Indicator for Use With Ethylene Oxide Sterilizers

Catalog No. EZEO

- Self-contained biological indicator for monitoring ethylene oxide sterilization

Color change is from red-to-yellow with growth

- Excellent indicator stability before and after exposure

48 hour results

- Process indicator label

ANSI/AAMI ST 79, class 1

ISO 11140-1, class 1

- Biological indicator

ANSI/AAMI ISO 11138

ORDERING INFORMATION

PART #	PACKAGED
EZEO	BOX = 100
EZEO/25	BOX = 25