

AUTOCLAVE TAPE



Steam/Chemical Vapor/Blue/EO Tapes

Available in Various Sizes

Catalog No. CI 122, CI 123, CI 124

- Sterilization package sealing tapes with a process indicator printed on the tape
 - Steam Tape indicator turns from white-to-brown/black when exposed to steam sterilization conditions. EO Gas process indicator turns from brown-to-green when exposed to ethylene oxide sterilization conditions
- Tape adheres to item throughout the sterilization process
- Shelf life:
 - Steam - 3 years from manufacture date
 - EO - 2 years from manufacture date.
- ANSI/AAMI ST79, class 1. ANSI/AAMI/ISO 11140-1, class 1

ORDERING INFORMATION		
PART #	DESCRIPTION	PACKAGED
CI 122/.5	1/2" STEAM TAPE	1 CASE = 72 ROLLS
CI 122/.75	3/4" STEAM TAPE	1 CASE = 48 ROLLS
CI 122/1	1" STEAM TAPE	1 CASE = 36 ROLLS
CI 123/1	1" BLUE STEAM TAPE	1 CASE = 36 ROLLS
CI 124/.75	3/4" EO TAPE	1 CASE = 48 ROLLS
CI 124/1	1" EO TAPE	1 CASE = 36 ROLLS