

WASH-CHECKS[®] PRO



Wash-Checks[®] Pro Rapid Test for Detecting Residual Protein

Catalog No. WC 112

- Wash-Checks Pro has the ability to detect residual proteins left behind on surfaces, endoscopes, and surgical instruments
 - Verifies the cleaning process
 - Versatile uses
- Can detect protein residues within one microgram sensitivity
 - If protein is present, the test will turn blue within ten seconds
 - Fastest results available
 - The higher the level of contamination of protein, the darker and faster the color will change to blue

ORDERING INFORMATION

PART #	PACKAGED
WC 112	1 CASE = 2 KITS OF 25 TESTS EACH



Please [click here](#) or scan the tag using the free mobile app from getscanlife.com to see our instructional video on the Wash-Checks Pro

SteriTec
www.steritecproducts.com