

WASH-CHECKS[®] PRO

Interpretation Guide

Wash-Checks[®] Pro Rapid Test

Catalog No. WC 112


UNEXPOSED

Protein Test liquid preliminary color is amber. The test will change to various shades of blue depending on the amount of protein detected.



HIGH LEVEL OF PROTEIN DETECTED

Protein Test liquid will change to various shades of blue depending on the amount of protein detected. Any shade of blue indicates that protein has been detected.



Unexposed

1µg

2µg

4µg

8µg

10µg

20µg

Positive Control