



Product Technical Data Sheet

Form No.: FTP-004 Rev.: G

Product Name: SteriTec Cross-Checks Indicator

Product Type: Chemical Indicator

Reorder Number: CI 107

Description: The SteriTec Cross-Checks is multi-variable indicator. It gives visual confirmation of the attainment of sterilization conditions in sterilizers operating at 132° – 134°C for 3 minutes or longer at the location of the indicator strip within the sterilizer chamber.

Performance Characteristics:

When all steam sterilization conditions are met, the check marks change color.

Color Change: White to Black

Critical Parameters: *(in a standard hospital steam sterilizer)*

3 minutes +/- .5 min at 132° - 134°C

Stated Values: *(as demonstrated in a steam resistometer)*

3 minutes (+0/-1.25 min) at 132°C (+0/-2°C)

2.5 minutes (+0/1.125 min) at 134°C (+0/-2°C)

Indicator Stability: Before and After Sterilization = Excellent

Physical Characteristics:

Dimensions: 0.5625" W x 7.66" L

Substrate Material: 80 pound Cover

Quantity / Bag: 250

Bags / Case: 16

Selling Unit: Case

Environmental Limitations:

Storage: 50° - 100°F (10° – 38°C) at 10 - 90% humidity

Shelf Life: The expiry date is printed on the product packaging.

FDA 510(k) Release to Market Date: August 19, 1996, K951113

Standards Classifications:

ANSI/AAMI ST79, Class 4

ISO 11140-1, Class 4