



Product Technical Data Sheet

Form No.: FTP-301 Rev.: C

Product Name: SteriTec Wash-Checks U
Cleaning Monitor

Product Type: Ultrasonic Washing Efficiency Monitor

SteriTec Reorder Number: WC 108, WC 108/R and WC 108/B

Description: The Wash-Checks U Monitor is an easy to use prepared test for determining the cleaning effectiveness of an ultrasonic cleaner. The test is placed in a holder (WC 102) for easy placement in the ultrasonic cleaner.

Performance Characteristics:

The test soil (blue square) simulates human blood and tissue debris. The removal of all test soil indicates effective cleaning.

Physical Characteristics:

The Wash-Checks U Monitor consists of a metal foil strip with a test soil applied (blue square).

Dimensions: 3.5" X 1"

Substrate Material: Aluminum Roll Stock

Units / Bag: 50

Bags / Case: 5

Selling Unit: Case

Environmental Limitations:

Storage: 50° - 100°F (10° - 38°C). Keep away from disinfectants and fumes.

Shelf Life: The expiry date is printed on the product packaging.

SteriTec FDA 510(k) Release to Market Date: N/A

Standards Classification: N/A