



Product Technical Data Sheet

Form No.: FTP-350 Rev.: B

Product Name: SteriTec Wash-Checks Pro Protein Detection Kit

Product Type: Washer Efficiency Monitor

SteriTec Reorder Number: WC 112, WC 112/R & WC 112/B

Description: The Wash-Checks Pro is a rapid test that has the ability to detect residual proteins left behind on the surfaces of washer-disinfectors, ultrasonic cleaners, endoscopes and surgical instruments.

Performance Characteristics:

This protein test can detect protein residues within 1 μg sensitivity and 10 seconds of time. When protein is present, the solution will turn from brown to blue. This test provides a semi-quantitative measure of surface cleanliness; the higher the level of contamination of protein (bioburden), the darker and faster the color changes to blue.

Physical Characteristics:

The SteriTec Wash-Checks Pro Protein Detection Kit contains 25 vials of protein detection solution, 2 Wash-Checks H monitors and 25 swabs.

Units / Box: 25

Boxes/Case: 2

WC 112/R : 8 Boxes/Case

Environmental Limitations:

Storage: Keep away from light and heat. Detection kit can be stored for up to 6 months at room temperature, or for up to 24 months under refrigeration.

Shelf Life: Shelf life is 2 years from manufacture of reagent when refrigerated and 6 months when not refrigerated.

SteriTec FDA 510(k) Release to Market Date: N/A

Standards Classification: N/A